

## Product Datasheet

# NATIVE Human APOLIPOPROTEIN A1 Protein (orb2022299)

<b>Description</b>	NATIVE Human APOLIPOPROTEIN A1 Protein
<b>Tested Applications</b>	ELISA
<b>Preservatives</b>	Lyophilized. 0.05 M Sodium chloride, 0.01 M Sodium carbonate
<b>Form/Appearance</b>	Lyophilized. 0.05 M Sodium chloride, 0.01 M Sodium carbonate
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Application notes</b>	Application Info: ntigen can be used for coating microplates and as a ligand for immunosorbent preparation.
<b>Purity</b>	SDS PAGE: > 90%
<b>MW</b>	28 kD
<b>Uniprot ID</b>	<a href="#">P02647</a>
<b>NCBI</b>	<a href="#">NP_000030.1</a>
<b>Expiration Date</b>	6 months from date of receipt.