

Product Datasheet

DHX9 antibody (orb20223)

Description	Goat polyclonal antibody to DHX9
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	DHX9 / RHA
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx. 140kDa band observed in nuclear lysates of cell line HeLa (calculated MW of 141kDa according to NP_001348.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	TEGRNALIHKSSVNC
Clonality	Polyclonal
MW	141
Entrez	1660
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20223
(1µg/ml)
staining of
nuclear He...



orb20223 (1
µg/mL)
staining of
nuclear ...