

Product Datasheet

NATIVE Rabbit ACE Protein (orb2022303)

| | |
|----------------------------|---|
| Description | NATIVE Rabbit ACE Protein |
| Tested Applications | FA |
| Preservatives | Purified protein from rabbit lung - liquid |
| Form/Appearance | Purified protein from rabbit lung - liquid |
| Storage | Store at -20°C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. |
| Note | For research use only |
| MW | Reported as approximately 129 kDa (sedimentation equilibrium with and without 6 M guanidine hydrochloride) and approximately 140 kDa (SDS gel electrophoresis and by gel filtration). |
| Uniprot ID | P12822 |
| NCBI | NP_001075864.1 |
| Expiration Date | 6 months from date of receipt. |