

Product Datasheet

PSPH antibody (orb20224)

Description	Goat polyclonal antibody to PSPH
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PSPH
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 26kDa band observed in Human Brain (Hippocampus) and in Mouse brain lysates (calculated MW of 25.0kDa according to Human NP_004568.2 and 25.1kDa according to Mouse NP_598661.1). Recommended concentration: 0.05-1 µg/ml.
Protein Sequence	RQQVKDNAKWYITD
Clonality	Polyclonal
MW	25
Entrez	5723
Dilution Range	ELISA: 1:32000, WB: 0.05-1 µg/ml
Expiration Date	12 months from date of receipt.



orb20224
(0.05ug/ml)
staining of
Mouse B...