

Product Datasheet

Mouse IgG3 Kappa (Îº) isotype Control [Out of stock] (orb2022792)

Description	Mouse IgG3 Kappa (Îº) isotype Control [Out of stock]
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, FLISA, IF, WB
Preservatives	Liquid (sterile filtered)
Form/Appearance	Liquid (sterile filtered)
Storage	2°C to 8°C -20°C
Note	For research use only
Application notes	Application Info: Mouse IgG3 kappa isotype control can be utilized as a control or standard reagent in Flow cytometry, Western Blotting, and ELISA experiments where determination of sample isotype is important.
Clonality	Monoclonal
Expiration Date	6 months from date of receipt.