

Product Datasheet

H1N1 California Protein (orb2023663)

Description	H1N1 California Protein
Species/Host	Virus
Preservatives	Liquid 10mM Sodium phosphate (pH 7), 150mM Sodium Chloride and 0.005% Tween 20
Form/Appearance	Liquid 10mM Sodium phosphate (pH 7), 150mM Sodium Chloride and 0.005% Tween 20
Storage	The Recombinant H1N1 A/California/04/2009 Recombinant should be stored at 4°C. DO NOT FREEZE
Note	For research use only
Purity	Greater than 90.0% under the conditions that would preserve its biological activity and tertiary structure.
MW	72 kDa
Expiration Date	6 months from date of receipt.