



Product Datasheet

GJA1 antibody (orb20237)


Description	Goat polyclonal antibody to GJA1
Species/Host	Goat
Reactivity	Canine, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	Connexin 43 / GJA1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 40kDa band observed in Rat Brain lysates (calculated MW of 43.0kDa according to Human NP_000156.1 and to Rat NP_036699.1). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	QPFDPPDDNQNSKK
Clonality	Polyclonal
MW	43
Entrez	2697
Dilution Range	ELISA: 1:32000, WB: 0.1-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human and Rat H...



Western blot analysis of Human and Rat H...



orb20237 (0.1 µg/mL) staining of Human ...