

Product Datasheet

H9N2 Hong-Kong Protein [Out of stock] (orb2024436)

Description	H9N2 Hong-Kong Protein [Out of stock]
Species/Host	Virus
Preservatives	The Recombinant H9N2 A/Hong-Kong/1073/99 solution contains 10mM Sodium Phosphate, pH 7.0 and 150mM NaCl and 0.01% Tween-20. Physical appearance: Sterile Filtered colorless solution.
Form/Appearance	The Recombinant H9N2 A/Hong-Kong/1073/99 solution contains 10mM Sodium Phosphate, pH 7.0 and 150mM NaCl and 0.01% Tween-20. Physical appearance: Sterile Filtered colorless solution.
Storage	A/H9N2 Hong-Kong/1073/99 Recombinant should be stored at 4C.
Note	For research use only
Purity	Greater than 90.0% as determined by SDS-PAGE.
Uniprot ID	Q9ICY5
NCBI	NP_859037.1
Expiration Date	6 months from date of receipt.