

Product Datasheet

H5N1 Vietnam Protein (orb2024439)

Description	H5N1 Vietnam Protein
Species/Host	Virus
Preservatives	The Recombinant H5N1 A/Vietnam/1203/04 solution contains 10mM Sodium Phosphate, pH 7.1, 0.005% Tween-20 and 150mM NaCl. Physical appearance: Sterile Filtered colorless solution.
Form/Appearance	The Recombinant H5N1 A/Vietnam/1203/04 solution contains 10mM Sodium Phosphate, pH 7.1, 0.005% Tween-20 and 150mM NaCl. Physical appearance: Sterile Filtered colorless solution.
Storage	H5N1 A/Vietnam/1203/04 Recombinant should be stored at 4C.
Note	For research use only
Purity	Greater than 90.0% as determined by Analysis by SDS-PAGE, HA1 and HA2 bands are observed using SDS-PAGE under reducing conditions.
Expiration Date	6 months from date of receipt.