

Product Datasheet

H1N1 New Caledonia Protein (orb2024445)

Description	H1N1 New Caledonia Protein
Preservatives	Sterile Filtered colorless solution FormulationThe H1N1 A/New Caledonia/20/99 IVR solution contains STE, 0.1 % sodium azide (NaN ₃) and 0.005 % thimerosal.
Form/Appearance	Sterile Filtered colorless solution FormulationThe H1N1 A/New Caledonia/20/99 IVR solution contains STE, 0.1 % sodium azide (NaN ₃) and 0.005 % thimerosal.
Storage	A/New Caledonia/20/99 IVR although stable 4C for 4 weeks, should be stored desiccated below -1°C. Please prevent freeze-thaw cycles.
Note	For research use only
Purity	Greater than 90.0% as determined byAnalysis by SDS-PAGE.
Expiration Date	6 months from date of receipt.