

Product Datasheet

H1N1 New Caledonia Protein (orb2024447)

Description	H1N1 New Caledonia Protein
Species/Host	Virus
Preservatives	Sterile Filtered colorless solution (1,666ug/ml). Formulation: The Recombinant H1N1 A/New Caledonia/20/99 solution 10mM Sodium phosphate, pH 7.2, 150mM NaCl.
Form/Appearance	Sterile Filtered colorless solution (1,666ug/ml). Formulation: The Recombinant H1N1 A/New Caledonia/20/99 solution 10mM Sodium phosphate, pH 7.2, 150mM NaCl.
Storage	H1N1 A/New Caledonia/20/99Recombinant should be stored at 4C.
Note	For research use only
Purity	Greater than 90.0% as determined byAnalysis by SDS-PAGE, HA1 and HA2 bands are observed using SDS-PAGE under reducing conditions.
Uniprot ID	Q6WG00
Expiration Date	6 months from date of receipt.