

Product Datasheet

Gliadin Wheat Protein (orb2025127)

Description	Gliadin Wheat Protein
Species/Host	E. coli
Preservatives	Sterile Filtered clear solution.
Form/Appearance	Sterile Filtered clear solution.
Storage	Store at 4C if entire vial will be used within 2-4 weeks. Store, frozen at -2°C for longer periods of time. Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Application Info: Western blot with monoclonal anti-hexa-His-tag antibody.
Purity	Greater than 94% as determined by SDS-PAGE.
Expiration Date	6 months from date of receipt.