

Product Datasheet

SAE1 Recombinant Protein (Human) [Out of stock] (orb2026163)

Description	SAE1 Recombinant Protein (Human) [Out of stock]
Species/Host	E. coli
Preservatives	Liquid. In HEPES, pH 8.0, containing 110mM potassium acetate, 2mM magnesium acetate, 0.5mM EGTA and 1mM DTT.
Form/Appearance	Liquid. In HEPES, pH 8.0, containing 110mM potassium acetate, 2mM magnesium acetate, 0.5mM EGTA and 1mM DTT.
Storage	Store at -8°C
Note	For research use only
MW	38 kDa, 90 kDa
Uniprot ID	O95605 , Q9UBE0
NCBI	NP_001014188.1
Expiration Date	6 months from date of receipt.