

## Product Datasheet

# CASP10 Recombinant Protein (Human) [Out of stock] (orb2026276)

<b>Description</b>	CASP10 Recombinant Protein (Human) [Out of stock]
<b>Species/Host</b>	E. coli
<b>Preservatives</b>	Liquid. In 50mM HEPES, pH 7.4, 100mM NaCl, 0.1% CHAPS, 0.09% APO-10, 1mM EDTA, 10% glycerol, 10mM DTT.
<b>Form/Appearance</b>	Liquid. In 50mM HEPES, pH 7.4, 100mM NaCl, 0.1% CHAPS, 0.09% APO-10, 1mM EDTA, 10% glycerol, 10mM DTT.
<b>Storage</b>	Store at -8°C
<b>Note</b>	For research use only
<b>Purity</b>	> 95% (SDS-PAGE)
<b>MW</b>	17 kDa, 12 kDa
<b>Uniprot ID</b>	<a href="#">Q92851</a>
<b>NCBI</b>	<a href="#">NP_001780.2</a>
<b>Expiration Date</b>	6 months from date of receipt.