

Product Datasheet

TG Native Protein (orb2026344)

Description	TG Native Protein
Conjugation	Unconjugated
Tested Applications	Control, WB
Preservatives	Freeze-dried Powder. PBS, pH 7.4, containing 0.01% SKL, 1mM DTT, 5% Trehalose and Proclin300.
Form/Appearance	Freeze-dried Powder. PBS, pH 7.4, containing 0.01% SKL, 1mM DTT, 5% Trehalose and Proclin300.
Storage	2°C to 8°C, -80°C
Tag	N-terminal His Tag
Note	For research use only
Purity	> 97%
Source	Native Extract
Uniprot ID	D1KKB3
NCBI	NP_001161890.1
Expiration Date	6 months from date of receipt.