

Product Datasheet

TPP1 antibody (orb20271)

Description	Goat polyclonal antibody to TPP1
Species/Host	Goat
Reactivity	Bovine, Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CLN2 / TPP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 65kDa band observed in Human Placenta lysates (calculated MW of 61.2kDa according to NP_000382.3). Recommended concentration: 0.3-1 µg/ml.
Protein Sequence	TPSVIRKRYNLTSQD
Clonality	Polyclonal
MW	61.2
Entrez	1200
Dilution Range	ELISA: 1:4000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Placenta ...