

## Product Datasheet

### TNFRSF1B Recombinant Protein (orb2027865)

<b>Description</b>	TNFRSF1B Recombinant Protein
<b>Species/Host</b>	E. coli
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Control, WB
<b>Preservatives</b>	Freeze-dried Powder. PBS, pH 7.4, containing 0.01% SKL, 5% Trehalose.
<b>Form/Appearance</b>	Freeze-dried Powder. PBS, pH 7.4, containing 0.01% SKL, 5% Trehalose.
<b>Storage</b>	2°C to 8°C -80°C
<b>Tag</b>	N-terminal His Tag
<b>Note</b>	For research use only
<b>Purity</b>	> 95%
<b>Source</b>	E.coli
<b>MW</b>	32kDa
<b>Uniprot ID</b>	<a href="#">P25119</a>
<b>NCBI</b>	<a href="#">NP_035740.2</a>
<b>Expiration Date</b>	6 months from date of receipt.