

Product Datasheet

TNFRSF1B Recombinant Protein (orb2027866)

Description	TNFRSF1B Recombinant Protein
Species/Host	E. coli
Conjugation	Unconjugated
Tested Applications	Control, WB
Preservatives	Freeze-dried Powder. PBS, pH 7.4, containing 0.01% SKL, 5% Trehalose.
Form/Appearance	Freeze-dried Powder. PBS, pH 7.4, containing 0.01% SKL, 5% Trehalose.
Storage	2°C to 8°C -80°C
Tag	N-terminal His Tag
Note	For research use only
Purity	> 95%
Source	E.coli
MW	32kDa
Uniprot ID	P25119
NCBI	NP_035740.2
Expiration Date	6 months from date of receipt.