



Product Datasheet

CXCL2 antibody (orb20284)

Description	Goat polyclonal antibody to CXCL2
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CXCL2 / MIP-2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 17Da band observed in Human Duodenum lysates (calculated MW of 11.4kDa according to NP_002080.1). Recommended concentration: 0.2-0.6 µg/ml.
Protein Sequence	KKIIEKMLKNGKSN
Clonality	Polyclonal
MW	11.4
Entrez	2920
Dilution Range	ELISA: 1:32000, WB: 0.03-0.1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Duodenum ...



Western blot analysis of Daudi lysate us...