

Product Datasheet

INDOL1 antibody (orb20285)

Description	Goat polyclonal to INDOL1
Species/Host	Goat
Reactivity	Bovine, Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	INDOL1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 50kDa band observed in Human Kidney lysates (calculated MW of 47.07kDa according to NP_919270.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	RDKTLESILHPR
Clonality	Polyclonal
MW	47.07
Entrez	169355
Dilution Range	ELISA: 1:64000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb20285
(0.5ug/ml)
staining of
Human Ki...