

Product Datasheet

CGA Recombinant Protein (orb2040609)

Description	CGA Recombinant Protein
Species/Host	E. coli
Conjugation	Unconjugated
Tested Applications	Control, WB
Preservatives	Freeze-dried Powder. PBS, pH 7.4, containing 0.01% SKL, 5% Trehalose.
Form/Appearance	Freeze-dried Powder. PBS, pH 7.4, containing 0.01% SKL, 5% Trehalose.
Storage	2°C to 8°C -80°C
Tag	N-terminal His Tag
Note	For research use only
Purity	> 90%
Source	E.coli
MW	14kDa
Uniprot ID	P01216
NCBI	NP_034019.1
Expiration Date	6 months from date of receipt.