


Product Datasheet

ZYX antibody (orb20416)

Description	Goat polyclonal antibody to ZYX
Species/Host	Goat
Reactivity	Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ZYX
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:4000.WB: Approx 85kDa band observed in Human Peripheral Blood Mononucleocytes lysates (calculated MW of 61.3kDa according to NP_003452.1). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources. Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	QQREKPRVQEKQH
Clonality	Polyclonal
MW	61.3
Entrez	7791
Dilution Range	ELISA: 1:4000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb20416
(0.1µg/ml)
staining of
Peripher...