


Product Datasheet

PDE1A antibody (orb20469)

Description	Goat polyclonal antibody to PDE1A
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	PDE1A
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 55kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 62.3kDa according to NP_005010.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	DLHKNSEDLVNAE
Clonality	Polyclonal
MW	62.3; 61.3
Entrez	5136
Dilution Range	ELISA: 1:64000, WB: 1-3 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



orb20469
(1ug/ml)
staining of
Human
Cere...