

Product Datasheet

APOB Recombinant Protein (orb2050092)

Description	APOB Recombinant Protein
Species/Host	E. coli
Conjugation	Unconjugated
Tested Applications	Control, WB
Preservatives	Freeze-dried Powder. PBS, pH 7.4, containing 0.01% SKL, 5% Trehalose.
Form/Appearance	Freeze-dried Powder. PBS, pH 7.4, containing 0.01% SKL, 5% Trehalose.
Storage	2°C to 8°C -80°C
Tag	N-terminal His Tag
Note	For research use only
Purity	> 95%
Source	E.coli
MW	24kDa
Uniprot ID	Q7TMA5
NCBI	NP_062160.2
Expiration Date	6 months from date of receipt.