

## Product Datasheet

### CLCN2 antibody (orb20502)

<b>Description</b>	Goat polyclonal antibody to CLCN2
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Bovine, Canine, Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	CLCN2
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 100kDa band observed in lysates of cell line HeLa (calculated MW of 98.5kDa according to NP_004357.3). Recommended concentration: 0.3-1 µg/ml.
<b>Protein Sequence</b>	QQLDEPVNFS DCK
<b>Clonality</b>	Polyclonal
<b>MW</b>	98.5
<b>Entrez</b>	<a href="#">1181</a>
<b>Dilution Range</b>	ELISA: 1:64000, WB: 0.3-1 µg/ml
<b>Expiration Date</b>	12 months from date of receipt.



orb20502  
(0.3ug/ml)  
staining of  
HeLa lys...