

Product Datasheet

CD3G Recombinant Protein (Cynomolgus monkey) [Out of stock] (orb2052775)

Description	CD3G Recombinant Protein (Cynomolgus monkey) [Out of stock]
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 6xHis-tagged
Note	For research use only
Protein Sequence	QSFEENRKLNVYNQEDGSVLLTCHVKNTNITWFKEGKMIDILTAHKNKWNLGSN TKDPRGVYQCKGSKDKSKTLQVYYRMCQNCIELNAAT
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	12.5 kDa
Uniprot ID	Q95LI7
Expiration Date	6 months from date of receipt.