

Product Datasheet

Alpha/beta-gliadin MM1 Recombinant Protein (Wheat) (orb2053226)

Description	Alpha/beta-gliadin MM1 Recombinant Protein (Wheat)
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Protein Sequence	VRVPVPLQLPQNPSQQQPQEQVPLVQQQFPGQQQFPFPQQYPQPQPFPSQ QPYLQLQPFPPQQLPYPQPQLPYPQPQLPYPQPQPFPPQQYPQSQPQYSQPQQ PISQQQQQQQQQQKQQQQQQQILQQILQQQLIPCRDVLQHQHSIAYGSS QVLQQSTYQLVQQLCCQQLWQIPEQSRCQAIHNVVHAILHQQQQQQQQQQ QPLSQVSFQQPQQQYPSGQGSFQPSQQNPQAQGSVQPQQLPFEEIRNLALET LPAMCNVYIPPYCTIAPVGIFGTN
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	49.3 kDa
Uniprot ID	P18573
Expiration Date	6 months from date of receipt.