

## Product Datasheet

# Alpha/beta-gliadin MM1 Recombinant Protein (Wheat) (orb2053227)

<b>Description</b>	Alpha/beta-gliadin MM1 Recombinant Protein (Wheat)
<b>Preservatives</b>	Liquid or Lyophilized powder
<b>Form/Appearance</b>	Liquid or Lyophilized powder
<b>Storage</b>	-20°C or -80°C
<b>Tag</b>	N-terminal 6xHis-SUMO-tagged
<b>Note</b>	For research use only
<b>Protein Sequence</b>	VRVPVPLQLPQNPSQQQPQEQVPLVQQQFPGQQQFPFPQQYPQPQPFPSQ QPYLQLQFPFPQQLPYPQPQLPYPQPQLPYPQPQPFPRQQYPYQSQPQYSQPQQ PISQQQQQQQQQQKQQQQQQQILQQILQQQLIPCRDVVLQQHSIAYGSS QVLQQSTYQLVQQLCCQQLWQIPEQSRCQAIHNVVHAILHQQQQQQQQQQQ QPLSQVSFQQPQQQYPSGQGSFQPSQQNPQAQGSVQPQQLPFEEIRNLALET LPAMCNVYIPPYCTIAPVGIFGTN
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Source</b>	E.coli
<b>MW</b>	49.3 kDa
<b>Uniprot ID</b>	<b>P18573</b>
<b>Expiration Date</b>	6 months from date of receipt.