

Product Datasheet

RBCL Recombinant Protein (Soybean) (orb2053232)

Description	RBCL Recombinant Protein (Soybean)
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Protein Sequence	PQTETKASVGFKAGVKDYKLYTTPDYETKDTDILAAFRVTPQPGVPPEEAGAAV AAESSTGTWTTVWTDGLTSLDRYKGRCYGLEPVAGEENQYIAYVAYPLDLFEEG SVTNMFTSIVGNVFGFKALRALRLEDLRIPTAYIKTFQGPPIHQVERDKLNKYGR PLLGCTIKPKLGLSAKNYGRAVYECLRGGLDFTKDDENVNSQPFMRWRDRFLFC AEAIFKSQAETGEIKGHYLNATAGTCEEMMKRAVFARELGVPIVMHDYLTGGFTA NTSLAHYCRDNGLLLIHHRAMHAVIDRQKNHGMHFRVLAKALRLSGGDHVVHAG TVVGKLEGEREITLGFVDLLRDDFVEKDRSRGIYFTQDWVSLPGVLPVASGGIHV WHMPALTEIFGDDSVLQFGGGTLGHPWGNAPGAVANRVALEACVQARNEGRD LAREGNEIIREASKWSPELAAACEVWKEIKFEFEAMDTL
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	68.4 kDa
Uniprot ID	P27066
NCBI	YP_538747
Expiration Date	6 months from date of receipt.