

Product Datasheet

VACWR062 Recombinant Protein (Vaccinia virus) (orb2053320)

Description	VACWR062 Recombinant Protein (Vaccinia virus)
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Protein Sequence	MDFIRRKYLIVTENNIDFLKDDTL SKVNNFTLNHVLALKYLVSNFPQH VITKDVLA NTNFFVFIHMRCKKVYEAVLRHAFDAP TLYVKALTKNYLSFSNTIQSYKETVHKL TQDEKFLEVAKYMDDELGELIGVNYDLVLNPLFHGG EPIKDM EIFFLKLFKKTDFKVV KKLSVIRLLIWAYLSKKDTGIEFADNDRQDIYTLFQQTGRIVHSNL TETFRDYIFPG DKTSYWVWLNESIANDADIVLNRHAITMYDKILSYIYSEIKQGRVNKNMLKLVYIF EPEKDIRELLEIIYDIPGDILSIIDAKNDDWK KYFISFYKANFINGNTFISDRTFNED LFRVVVQIDPEYFDNERIMSLFSTSAADIKRFDELDINNSYISNIIYEVNDITLDTMD DMKKCQIFNEDTSYYVKEYNTYLFLHESDPMVIENGILKKLSSIKSKSRRLNLFSKN ILKYYLDGQLARLGLVLD DYKGDLLVKMINHLKSVEDVSAFVRFSTDKNPSILPSLI KTILASYNISIVLFQRFLRDNLYHVEEFLDKSIHLTKTDKKYILQLIRHGRS
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	82.7 kDa
Uniprot ID	P21607
NCBI	YP_232944
Expiration Date	6 months from date of receipt.