

Product Datasheet

HA Recombinant Protein [Out of stock] (orb2054708)

Description	HA Recombinant Protein [Out of stock]
Species/Host	Virus
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Protein Sequence	DTICIGYHANNSTDTVDTVLEKNVTVTHSVNLLSDSHNGKLCRLKGIAPLQLGKC NIAGWILGNPECESLLSKRSWSYIAETPNSGTCYPGDFADYEELREQLSSVSS FERFEIFPKERSWPNNINIGVTAACSHAGKSSFYKNLLWLTEKDGSYPNLNKS VVKKEKEVLVLWGVHHPNINENQKTLRKENAYVSVSSNYNRRFTPEIAERP RGQAGRMNYYWTLLEPGDTIIFEANGNLIAPWYAFALSRGLGSGIITSNASMDEC DTKCQTPQGAINSSLPFQNIHPFTIGECPKYVRSTKLRMVTGLRNIPSIQS
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	52.6 kDa
Uniprot ID	Q0HD60
Expiration Date	6 months from date of receipt.