

## Product Datasheet

### X Recombinant Protein (orb2054787)

<b>Description</b>	X Recombinant Protein
<b>Species/Host</b>	Virus
<b>Preservatives</b>	Liquid or Lyophilized powder
<b>Form/Appearance</b>	Liquid or Lyophilized powder
<b>Storage</b>	-20°C or -80°C
<b>Tag</b>	N-terminal 6xHis-SUMO-tagged
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MAARLCCQLDPARDVLCRLPVGAE SRGRPFSGSLGTLSSPSPSAVPTDHGAHLS LRGLPVCAFSSAGPCALRFTSARRMETTVNAHQILPKVLHKRTLGLSAMSTTDLE AYFKDCLFKDWEELGEEIRLKVFLGGCRHKLVCAPAPCNFF TSA
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Source</b>	E.coli
<b>MW</b>	32.6 kDa
<b>Uniprot ID</b>	<b>P03165</b>
<b>Expiration Date</b>	6 months from date of receipt.