

## Product Datasheet

# C3L Recombinant Protein [Out of stock] (orb2055521)

<b>Description</b>	C3L Recombinant Protein [Out of stock]
<b>Species/Host</b>	Virus
<b>Preservatives</b>	Liquid or Lyophilized powder
<b>Form/Appearance</b>	Liquid or Lyophilized powder
<b>Storage</b>	-20°C or -80°C
<b>Tag</b>	N-terminal 6xHis-tagged
<b>Note</b>	For research use only
<b>Protein Sequence</b>	CCTIPSRPINMKFKNSVETDANANYNIGDTIEYLCLPGYRKQKMGPIYAKCTGTG WTLFNQCIKRRCPSPRIDIDNGQLDIGGVDFGSSITYSCNSGYHLIGESKSYCELG STGSMVWNPEAPICESVKCQSPPSISNGRHNGYEDFYTDGSSVVTYSCNSGYSLI GNSGVLCSSGGEWSDPPTCQIVKCPHTISNGYLSSGFKRSYSYNDNVDFKCKYG YKLSGSSSSTCSPGNTWKPELPCVCR
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Source</b>	Yeast
<b>MW</b>	28.6 kDa
<b>Uniprot ID</b>	<a href="#">P68638</a>
<b>NCBI</b>	<a href="#">YP_232907.1</a>
<b>Expiration Date</b>	6 months from date of receipt.