

Product Datasheet

UL111A Recombinant Protein (orb2055605)

| | |
|-------------------------|--|
| Description | UL111A Recombinant Protein |
| Species/Host | Virus |
| Preservatives | Liquid or Lyophilized powder |
| Form/Appearance | Liquid or Lyophilized powder |
| Storage | -20°C or -80°C |
| Tag | N-terminal 6xHis-SUMO-tagged |
| Note | For research use only |
| Protein Sequence | SEEAKPATTTIKNTKPQCRPEDYATRLQDLRVTFHRVKPTLQREDDYSVWLDGT VVKGCWGCVMDWLLRRYLEIVFPAGDHVYPGLKTELHSMRSTLESYKDMRQ CPLLGCGDKSVISRLSQEAERKSDNGTRKGLSELDTLFSRLEEYLHSRK |
| Purity | Greater than 90% as determined by SDS-PAGE. |
| Source | E.coli |
| MW | 34 kda |
| Uniprot ID | P17150 |
| Expiration Date | 6 months from date of receipt. |