

Product Datasheet

UL111A Recombinant Protein (orb2055607)

Description	UL111A Recombinant Protein
Species/Host	Virus
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 6xHis-SUMO-tagged
Note	For research use only
Protein Sequence	SEEAKPATTTIKNTKPQCRPEDYATRLQDLRVTFHRVKPTLQREDDYSVWLDGT VVKGCWGCVMDWLLRRYLEIVFPAGDHVYPGLKTELHSMRSTLESYKDMRQ CPLLGCGDKSVISRLSQEAERKSDNGTRKGLSELDTLFSRLEEYLHSRK
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	34 kda
Uniprot ID	P17150
Expiration Date	6 months from date of receipt.