


Product Datasheet

SLC7A7 antibody (orb2056)

Description	Goat polyclonal antibody to SLC7A7
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	SLC7A7 / γ+LAT-1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 55-58kDa band observed in Human Lung and Skeletal Muscle lysates (calculated MW of 56.0kDa according to NP_001119577.1. Recommended concentration: 1-3µg/ml.
Protein Sequence	SQPEVETSPLGD
Clonality	Polyclonal
MW	56
Entrez	9056
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Lung (Lan...