

Product Datasheet

ILDR2 Recombinant Protein (orb2056082)

Description	ILDR2 Recombinant Protein
Species/Host	Human
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 6xHis-tagged
Note	For research use only
Protein Sequence	MDRVLLRWISLFWLTAMVEGLQVTVPDKKKVAMLFQPTVLRCHFSTSSHQPAV VQWKFKSYCQDRMGESLGMSSTRAQSLSKRNLEWDPYLDCLDSRRTVRVVAS KQGSTVTLGDFYRGREITIVHDADLQIGKLMWGDGSLYYCIITTPDDLEGKNEDS VELLVLGRTGLLADLLPSFAVEIMPE
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	23.1 kDa
Uniprot ID	Q71H61
NCBI	NP_955383.1
Expiration Date	6 months from date of receipt.