

Product Datasheet

GP Recombinant Protein [Out of stock] (orb2056091)

Description	GP Recombinant Protein [Out of stock]
Species/Host	Virus
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 6xHis-tagged
Note	For research use only
Protein Sequence	QTNTKATGKCNPNLHYWTAQEQHNAAGIAWIPYFGPGAEGIYTEGLMHNQNAL VCGLRQLANETTQALQLFLRATTELRTYTILNRKAIDFLLRRWGGTCRILGPDCCI EPHDWTKNITDKINQIHDFFIDNPLPN
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	17.4 kDa
Uniprot ID	Q7T9D9
NCBI	YP_138523.1
Expiration Date	6 months from date of receipt.