

Product Datasheet

E6 Recombinant Protein [Out of stock] (orb2056279)

Description	E6 Recombinant Protein [Out of stock]
Species/Host	Virus
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 6xHis-tagged
Note	For research use only
Protein Sequence	MHQKRTAMFQDPQERPRKLPQLCTELQTTIHDIILECVYCKQQLLRREVDFAFR DLCIVYRDGNPYAVCDKCLKFYKISEYRHYCYSLYGTTLEQQYNKPLCDLLIRCIN CQKPLCPEEKQRHLDKKQRFHNIIRGRWTGRCMSSCRSSRTRRETQL
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	21.2 kDa
Uniprot ID	P03126
NCBI	NP_041325.1
Expiration Date	6 months from date of receipt.