

Product Datasheet

VP2 Recombinant Protein [Out of stock] (orb2056326)

Description	VP2 Recombinant Protein [Out of stock]
Species/Host	Virus
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 6xHis-tagged
Note	For research use only
Protein Sequence	GAALALLGDLVASVSEAAAATGFSVAEIAAGEAAAAIEVQIASLATVEGITSTSEAI AAIGLTPQTYAVIAGAPGAIAGFAALIQTVSGISSLAQVGYRFFSDWDHKVSTVGL YQQSGMALELFNPDEYYDILFPGVNTFVNNIQYLDPRHWGPSLFATISQALWHVI RDDIPSITSQELQRRTERFFRDSLARFLEETTWTIVNAPINFYNIQQYSDLSPIRP SMVRQVAEREGTRVHFGHTYSIDDADSIEEVTQRMDLRNQQSVHSGEFIEKTIA PGGANQRTAPQWMLPLLLGLYGTVTPALEAYEDGPNQKRRVSRGSSQKAKGT RASAKTTNKRRSRSSRS
Purity	Greater than 90% as determined by SDS-PAGE.
Source	Yeast
MW	40.2 kDa
Uniprot ID	P03094
NCBI	YP_717937.1
Expiration Date	6 months from date of receipt.