


Product Datasheet

NDUFS2 antibody (orb20587)

Description	Goat polyclonal antibody to NDUFS2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NDUFS2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 48kDa band observed in Human, Mouse and Rat Heart lysates (calculated MW of 51.9kDa according to NP_001159631.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	SRGVRQWQPDVE
Clonality	Polyclonal
MW	52.5; 51.9
Entrez	4720
Dilution Range	ELISA: 1:16000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart lys...