



Product Datasheet

PTPRB antibody (orb20615)

Description	Goat polyclonal to PTPRB
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	PTPRB
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms (NP_001103224.1; NP_002828.3).ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Preliminary experiments in Human, Mouse and Rat Testis lysates gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	ENHSFSQERTVDPK
Clonality	Polyclonal
MW	249; 224; 214; 214
Entrez	5787
Dilution Range	ELISA: 1:16000, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A431 cell...



Flow cytometric analysis of A431 cells u...