

Product Datasheet

L2RB Recombinant Protein (Mouse) [Out of stock] (orb2062375)

Description	L2RB Recombinant Protein (Mouse) [Out of stock]
Preservatives	Liquid or Lyophilized powder
Form/Appearance	Liquid or Lyophilized powder
Storage	-20°C or -80°C
Tag	N-terminal 6xHis-tagged
Note	For research use only
Protein Sequence	ASAAVKNCSHLECFYNSRANVSCMWSHEEALNVTTCHVHAKSNLRHWNKTCEL TLVRQASWACNLILGSFPESQSLTSVDLLDINVVCWEEKGWRRVKTCDFHFPDN LRLVAPHSLQVLHIDTQRCNISWKVSQVSHYIEPYLEFEARRLLGHSWEDASVL SLKQRQQWLFLEMLIPSTSYEVQVRVKAQRNNTGTWSPWSQPLTFRTRPADPM KE
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	29.3 kDa
Uniprot ID	P16297
Expiration Date	6 months from date of receipt.