


Product Datasheet

NUP107 antibody (orb20630)

Description	Goat polyclonal to NUP107
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	NUP107
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:8000.WB: Approx 100kDa band observed in Human Peripheral Blood Mononucleocytes lysates (calculated MW of 106kDa according to NP_065134.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	PLPAEDDNAIREH
Clonality	Polyclonal
MW	106
Entrez	57122
Dilution Range	ELISA: 1:8000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Peripheral Bloo...



Western blot analysis of nuclear lysate ...