

## Product Datasheet

# Rabbit anti-Sheep Cell Hemolysin [Out of stock] (orb2063905)

<b>Description</b>	Rabbit anti-Sheep Cell Hemolysin [Out of stock]
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	HAI
<b>Immunogen</b>	Sheep washed pooled Red Blood Cells (RBC)
<b>Preservatives</b>	Lyophilized
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Application Info: Rabbit Sheep Cell Hemolysin is suitable for agglutination, complement fixation related assays.
<b>Expiration Date</b>	12 months from date of receipt.