



Product Datasheet

HOXC8 antibody (orb20656)

Description	Goat polyclonal to HOXC8
Species/Host	Goat
Reactivity	Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	HOXC8
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 30kDa band observed in Human Peripheral Blood Mononucleocytes lysates (calculated MW of 27.8kDa according to NP_073149.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	SLYGAQQEAS
Clonality	Polyclonal
MW	27.8
Entrez	3224
Dilution Range	ELISA: 1:2000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Peripheral Blood Mononucleocytes lysates



Western blot analysis of Mouse Spinal Cord lysates