

Product Datasheet

ACOX3 antibody (orb20682)

Description	ACOX3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Peroxisomal acyl-coenzyme A oxidase 3 protein (401-624AA)
Target	Major pollen allergen Dac g 4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15254
Dilution Range	WB:1:500-5000, IHC-P:1:10-100, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.