

Product Datasheet

ACOX3 antibody (HRP) (orb20683)

Description	ACOX3 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Peroxisomal acyl-coenzyme A oxidase 3 protein (401-624AA)
Target	Major pollen allergen Dac g 4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15254
Dilution Range	WB:1:500-5000, IHC-P:1:1-10, IF:1:50-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.