

Product Datasheet

HSF1 Antibody (Phospho-Ser326) : FITC [Out of stock] (orb2069046)

Description

HSF1 Antibody (Phospho-Ser326) :

Species/Host	Rabbit
Conjugation	FITC
Tested Applications	IF, IHC-P, WB
Immunogen	KLH conjugated synthetic phosphopeptide derived from human HSF1 around the phosphorylation site of Ser326
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer.
Storage	All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.
Note	For research use only
Application notes	Application Info: WB(1:100-1000)IHC-P(1:100-500)IF(1:50-200)
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q00613
Expiration Date	12 months from date of receipt.