

---

## Product Datasheet

**TP53BP1 Antibody (Phospho-Ser25/Ser29) : Biotin  
[Out of stock] (orb2069257)**

**Description**

TP53BP1 Antibody (Phospho-

stock]

**Species/Host**

Rabbit

**Conjugation**

Biotin

**Tested Applications**

IF, IHC-P

**Immunogen**

KLH conjugated synthetic phosphopeptide derived from human p53BP1 around the phosphorylation site of Ser25/29

**Preservatives**

Liquid. Purified antibody supplied in 1x PBS buffer.

**Form/Appearance**

Liquid. Purified antibody supplied in 1x PBS buffer.

**Storage**

All conjugated antibodies should be stored in light-protected vials or covered with a light protecting material (i.e. aluminum foil). Conjugated antibodies are stable for at least 12 months at 4C. If longer storage is desired (24 months), conjugates may be diluted with up to 50% glycerol and stored at -2°C to -8°C. Freezing and thawing conjugated antibodies will compromise enzyme activity as well as antibody binding.

**Note**

For research use only

**Application notes**

Application Info: IHC-P(1:100-500)IF(1:50-200)

**Isotype**

IgG

**Clonality**

Polyclonal

**Uniprot ID**
[Q12888](#)
**Expiration Date**

12 months from date of receipt.